



【Product Name】

Tumour Treating Fields Instrument

【Product Introduction】

Tumour treating fields instrument applies an alternating electric field through insulated electrode patches to prevent the formation of spindle microtubules in tumour cells and the separation of intracellular organelles during cell division phase, and induce apoptosis in the division phase. It has the advantages of safety and little side effects, and has the potential of treating many kinds of malignant tumours. It can be applied to the treatment of a variety of solid tumours to prolong the progression free survival (PFS) and overall survival (OS) of tumour patients.