



【Product Name】

Drug-Eluting PTA Balloon Dilatation Catheters

【Product Introduction】

This product is an OTW balloon dilatation catheter, which consists of an terminal, a balloon development site, a balloon with drug coating, a catheter and a catheter holder. The drug coating is pure paclitaxel. The paclitaxel coating is evenly distributed on the surface of the effective length of the balloon, and the drug dosage is $3 \mu\text{g}/\text{mm}^2 \pm 1 \mu\text{g}/\text{mm}^2$. The product is sterilized with ethylene oxide and is disposable after using. It is suitable for percutaneous transluminal angioplasty in patients with stenosis or occlusion of the femoral and popliteal arteries (except infrapopliteal artery).