



【Product Name】

Silk Fibroin Wound Dressing (Sidaiyi)

【Product Introduction】

Sidaiyi is a disposable medical biological sterile wound dressing made from natural silk fibroin. It adopts a unique double-layer structure design, with the upper layer of silicone rubber membrane, which can block bacteria and maintain the wet environment of the wound, and the lower layer of porous silk fibroin membrane, which can guide the regeneration of tissues and absorb the wound exudate. The product has excellent flexibility and compliance. Under the wet state, it can be very good compliance with the wound tissue, and can be used for a variety of clinical wound tissue repair.