



【Product Name】

Mitral Valve Repair Device

【Product Introduction】

The product, based on surgical mitral valve edge-to-edge repair technology, utilizes mechanical clamping principle to implanted in the human body through the transcatheter implantation, then the upper clamp, lower clamp of the mitral valve clamping device clamp the mitral valve leaflets, form a "V"-shape with the two arms of the clip converging towards the midline, and at the same time, drive the clamped tissues to the midline. Once enclosed within the closure ring, the mitral valve is locked and then released. The implantation of the mitral valve clips can reduce mitral regurgitation and achieve the effect of treating regurgitation.