



**【Product Name】**

Transcatheter Aortic Valve Intervention

**【Product Introduction】**

Under the supervision of medical imaging equipment, the product is implanted through femoral artery via a catheter into the aortic annulus of the human body, replacing the original diseased aortic valve, to improve the stenosis of the diseased area, and enhance cardiac function. It is indicated for patients who after evaluation by a heart team using a scoring system, are deemed to have symptomatic, calcified, severe degenerative native aortic valve stenosis and are unsuitable for conventional surgical valve replacement, particularly those aged 70 years or older.