



【Product Name】

Annuloplasty Ring

【Product Introduction】

The products can restore mitral or tricuspid annulus to normal size and effectively promote the remodelling of the valve's physiological structure for a long term and better adaptation to the atrioventricular annular activity of the cardiac cycle. The saddle-shaped annulus design is used in most cases, and it can be classified into direct and indirect annuloplasty techniques based on the device and surgical approach used. It is suitable for patients who need mitral and tricuspid valve repair due to various atrioventricular valve diseases, such as mitral regurgitation, tricuspid regurgitation, etc.