



【Product Name】

Implantable Neurostimulation System

【Product Introduction】

The products are designed to stimulate target nerves with a certain level of current pulses in order to adjust or restore the function of brain, nerves or muscles, thus improving the symptoms and quality of life of diseased people. It usually consists of an implantable pulse generator and accessories. By applying electrical impulses to specific areas of brain or nervous system, it treats diseases such as Parkinson's disease, controls epilepsy, alleviates chronic intractable pain in the trunk and/or limbs, aids in intestinal and urinary control, and addresses dystonia, among others.